

# Enteral Feeding Threeway Stopcock



**Product Categories:** [Enteral Devices](#)  
**Product Page:** <https://www.medzus.com/product/enteral-feeding-threeway-stopcock/>

## Product Summary



Do Not Re-Sterilize



Sterilized Using Ethylene Oxide



Latex Free

## Product Description

- Neuraxial male & female lock connector (NRfit: ISO 80369-6) - A yellow colour rotator for ease of neuraxial colour identification **Available Configuration** - Also available Lipid Resistance & BPA free

## Product Attributes

- Brand Name: Neureo-Fit
- Product Code: NR-6003
- Symbols: Latex Free, Sterile EO, Do-Not resterilize.

## Product Gallery



