

# Enteral Feeding Bag



**Product Categories:** [Enteral Devices](#)

**Product Page:** <https://www.medzus.com/product/enteral-feeding-bag/>

## Product Summary



Latex Free



Sterilized Using Ethylene Oxide



Do Not Re-Sterilize

## Product Description

- Enteral male & female lock connectors (ENfit : ISO 80369-3) - A purple colour components for ease of enteral colour identification - Enteral feeding for Adult, Pediatric & Neonate with safe connection; ISO 80369-3 compliant connectors - Rear side of bag provided with pouch for ice to preserve the food for longer intervals of time - A hanger is provided with bag for easy handling - Smooth roller clamp facilitates easy, safe control and adjustment of flow rate. - Clear Drip Chamber to see the flow - ENfit male lock Y-connector for intermittent administration - Symbols- Latex Free, Sterile EO, Do-Not resterilize.

### Available Configuration-

- Feeding Bag 500/1000/1200 ml - Customization Available - - PVC Tubing (DEHP/Non-DEHP) with Customized Length - Silicon Tubing (Optional) - Silicon Connector (Optional)

## Product Attributes

- Brand Name: Entero-Fit
- Product Code: ENT-3006