

# Enteral Feeding Set



**Product Categories:** [Enteral Devices](#)

**Product Page:** <https://www.medzus.com/product/enteral-feeding-set/>

## Product Summary



Do Not Re-Sterilize



Sterilized Using Ethylene Oxide



Latex Free

## Product Description

Enteral Feeding Set Comprised of-

- Enteral Spike with Cap - PVC Tubing (DEHP/Non-DEHP) - Roller Clamp - Drip Chamber - Silicon Tubing (Optional) - Silicon Connector (Optional) - Feeding Y-Connector with ENfit Female Luer Lock with Cap - ENfit Female Luer Lock with Cap - Kink Resistance Tubing - A purple colour components for ease of enteral colour identification - Enteral male & female lock connectors (ENfit : ISO 80369-3) - Symbols- Latex Free, Sterile EO, Do-Not resterilize.

Available Configuration-

- PVC Tubing (DEHP/Non-DEHP) with Customized Length - Silicon Tubing (Optional) - Silicon Connector (Optional)

## Product Attributes

- Brand Name: Entero-Fit
- Product Code: ENT-3005