

Enteral Feeding Threeway Stopcock



Product Categories: [Enteral Devices](#)

Product Page: <https://www.medzus.com/product/enteral-feeding-threeway-stopcock/>

Product Summary



Do Not Re-Sterilize



Sterilized Using Ethylene Oxide



Latex Free

Product Description

- Neuraxial male & female lock connector (NRfit: ISO 80369-6) - A yellow colour rotator for ease of neuraxial colour identification **Available Configuration** - Also available Lipid Resistance & BPA free

Product Attributes

- Brand Name: Neureo-Fit
- Product Code: NR-6003
- Symbols: Latex Free, Sterile EO, Do-Not resterilize.

Product Gallery



