

# Enteral Feeding Tube



Product Categories: [Enteral Devices](#)

Product Page: <https://www.medzus.com/product/enteral-feeding-tube/>

## Product Summary



Do Not Re-Sterilize



Sterilized Using Ethylene Oxide



Latex Free

## Product Description

- Enteral male lock connector (ENfit : ISO 80369-3) - A purple colour male lock connector for ease of enteral colour identification - Size colour coded identification - For neonatal and paediatric nutritional feeding - DEHP free with graduation mark after every cm - Closed distal end with two lateral eyes - Soft, frosted & kink resistant PVC tubing with radio-opaque line - Radio opaque line to help confirm tube placement and migration monitoring

French Size  
Colour Code

FG-4  
Red

FG- 5  
Grey

FG- 6  
Light Green

FG-7  
White

FG- 8  
Blue

FG-9  
Orange

FG 10  
Black

## Product Attributes

- Brand Name: Entero-Fit
- Product Code: ENT-3003.8.40
- Symbols: Latex Free, Sterile EO, Do-Not resterilize.
- Color: Black, Blue, Grey, Light Green, Orange, Red, White
- Sizes: FG-10, FG-4, FG-5, FG-6, FG-7, FG-8, FG-9

## Product Gallery

