

Enteral Feeding Tube



Product Categories: [Enteral Devices](#)
Product Page: <https://www.medzus.com/product/enteral-feeding-tube/>

Product Summary



Do Not Re-Sterilize



Sterilized Using Ethylene Oxide



Latex Free

Product Description

- Enteral male lock connector (ENfit : ISO 80369-3) - A purple colour male lock connector for ease of enteral colour identification - Size colour coded identification - For neonatal and paediatric nutritional feeding - DEHP free with graduation mark after every cm - Closed distal end with two lateral eyes - Soft, frosted & kink resistant PVC tubing with radio-opaque line - Radio opaque line to help confirm tube placement and migration monitoring

French Size
Colour Code

FG-4
Red

FG- 5
Grey

FG- 6
Light Green

FG-7
White

FG- 8
Blue

FG-9
Orange

FG 10
Black

Product Attributes

- Brand Name: Entero-Fit
- Product Code: ENT-3003.8.40
- Symbols: Latex Free, Sterile EO, Do-Not reesterilize.
- Color: Black, Blue, Grey, Light Green, Orange, Red, White
- Sizes: FG-10, FG-4, FG-5, FG-6, FG-7, FG-8, FG-9

Product Gallery

