

Luer to ENfit Transition Adaptor



Product Categories: [Enteral Devices](#)
Product Page: <https://www.medzus.com/product/luer-to-enfit-transition-adaptor/>

Product Summary



Do Not Re-Sterilize



Sterilized Using Ethylene Oxide



Latex Free

Product Description

- Enteral female lock (ENfit : ISO 80369-3) with transition from female luer lock connector - A purple colour for ease of enteral colour identification

Product Attributes

- Brand Name: Entero-Fit
- Product Code: ENT-3007