

# Luer to ENfit Transition Adaptor



**Product Categories:** [Enteral Devices](#)  
**Product Page:** <https://www.medzus.com/product/luer-to-enfit-transition-adaptor/>

## Product Summary



Do Not Re-Sterilize



Sterilized Using Ethylene Oxide



Latex Free

## Product Description

- Enteral female lock (ENfit : ISO 80369-3) with transition from female luer lock connector - A purple colour for ease of enteral colour identification

## Product Attributes

- Brand Name: Entero-Fit
- Product Code: ENT-3007