

# Ommaya Reservoir Kit



NO IMAGE AVAILABLE

**Product Categories:** [Neuraxial Devices](#)

**Product Page:** <https://www.medzus.com/product/ommaya-reservoir-kit/>

## Product Description

- All components connection is compatible to neuraxial small bore connectors (NRfit: ISO 80369-6) - Kit Consist of - - Neuraxial Syringe - Scalp Vein Set - Safety needle shield minimizes the risk of needle stick injuries. - Color coded flexible wings for easy identification of size. - Butterfly shaped wings facilitate easy handling and attachment with the skin. Butterfly is connected to soft non toxic, non-irritant, medical grade PVC tubing. - The product is also available with DEHP free tubing.

## Product Attributes

- Brand Name: Neureo-Fit
- Product Code: NR-6007
- Sizes: 18G to 25G